

Table 1. List of primers and probes used for real-time RT-PCR.

GENE	PUBMED ID	Forward primer	Reverse primer	Probe
Bcl-2	M14745	CATGTGTGTGGAGAGCGTCAA	AACAGAGGCCGCATGCTG	CCTGGTGGACAACATCGCCCTGT
Bcl-xL	Z23115	CCACTTACCTGAATGACCACCTAGA	CTGGCCCTTTTCGGCTCTC	TGCTGCATTGTTCCCATAGAGTTCCACA
Bcl-w	U59747	GATGGTGGCCTACCTGGAGA	TGTCCTCACTGATGCCAGTT	ACTGCTGTGGATCCAGTCAGCCAGCCGC
Bax	L22473	AGTAACATGGAGCTGCAGAGGAT	GCTGCCACTCGGAAAAAGAC	TTGCCGCCGTGGACACAGACTC
Bcl-2A1	NM_004049	ACACAGGAGAATGGATAAAGGCAA	CATCCAGCCAGATTTAGGTTCAA	TTACAAAGCCATTTTCCAGCCTCCGT
MCL-1, variant 1	NM_021960	GCCTTTGTGGCTAAACACTTGA	ACAAACCCATCCCAGCCTCT	ACCATAAACCAAGAAAGCTGCATCGAACC
MCL-1, variant 2	NM_182763	GCAGCGCAACCACGAGA	ATGCCACCTTCTAGGTCCTCTACA	TGGAAGAACTCCACAAACCCATCCTTGG